OralScan AI

Agentic AI for Oral Cancer Screening & Tooth Health

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User Guide & Technical Report — Bilingual Edition

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Important: Educational prototype — NOT a medical device, NOT a diagnosis. Always consult a qualified dentist or doctor.

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ONE photo returns TWO results: Function 1 — oral-cancer cues (visual finding, region, probability, 0–100 risk score) and Function 2 — tooth health (cavity 蛀牙 risk, gum condition, plaque, staining, hygiene score, and whether professional scaling is needed).

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On Medium/High risk a 3-expert AI panel automatically debates the case and a moderator summarises the verdict. The History page shows your risk-trend chart.

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3. Doctor — Review & AI Agents • 医生端：审阅与 AI 智能体

- Triage Queue: every screening from web AND Telegram, sorted by urgency.
分诊队列：来自网页和 Telegram 的所有筛查，按紧急程度排序。

- Patients page: click ' Load demo patients' once — Aisyah / Tan / Muthu appear with full records (name, age, race, habits, brushing times, last scaling date, prior reports, charts, heatmap).

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- AI Triage Agent (top of Patients page): ask 'Which patient should I see first?' — the agent calls its compare_all_patients tool and ranks the cohort with reasons. Watch the tool-trace chips.

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- Per-patient AI Assistant: summarise history, risk trend, red flags, scaling status — answers with real numbers fetched live by its tools.

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- Session review: agree / upgrade / downgrade / refer, with notes — saved to the record.

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Log in as admin / 123456 → you land on /admin. Paste any of the 3 provider keys (ChatGPT / Gemini / Claude) and click Save. They apply to EVERY user instantly — web and the Telegram bot — no redeploy. Keys are masked after saving.

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Key precedence everywhere: user's own pasted key → admin key → backend env key → offline mock.

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5. Telegram Bot — Zero-Press • Telegram 机器人：零点击

Just send the bot a photo of the inside of a mouth — that's it. No buttons, no questions. The full agent pipeline runs immediately using the active demo patient's stored records (history, habits, last scaling date) and replies with BOTH reports plus the agent trace.

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6. Two-Minute Demo Script • 两分钟演示脚本

- ① Login doctor/doctor → Patients → 'Load demo patients' → ask the triage agent who to see first.

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- ② Login patient1 → New Screening → run a sample case → show the two-function result and the live agent pipeline animation.

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- ③ On Telegram, send a photo → instant result; /change patient → send again as high-risk Tan.

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- ④ Open /introduction → walk the architecture diagram and the OpenClaw / Claude Code / Hermes pattern table.

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PART 2 • TECHNICAL REPORT — 技术报告

1. System Overview • 系统概览

OralScan AI is a hierarchical multi-agent system: ~14 specialised agents collaborate on every screening. One photo enters through any channel (web or Telegram); a single gateway routes it to the Orchestrator, which dispatches the specialist agents in order; every LLM output is structured JSON; every result becomes a shared record visible to both patient and doctor.

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3. The Agent Registry • 智能体清单

Agent 智能体	Role 职责
Orchestrator 编排者	Knows the full pipeline; dispatches every agent; writes the audit log. 掌握全流程，逐一调度并记录审计日志。
Toothbrush IoT 牙刷物联网	Simulated smart-toothbrush telemetry: duration, pressure, coverage. 模拟智能牙刷遥测：时长、压力、覆盖率。
Vision Screening 视觉筛查	ONE image → TWO functions: cancer cues + tooth-health signs; provider pluggable with mock fallback. 一图双检：癌症线索 + 牙齿健康；模型可插拔并自动回退。
Cancer Risk Scoring 风险评分	Deterministic 0–100 score, literature-calibrated weights (IARC / INHANCE). 确定性 0–100 评分，权重源自文献。
Dental Wellness 牙齿健康	Cavity 蛀牙 / gums / plaque / staining / hygiene + scaling (洗牙) advice using the last-scaling record. 结合上次洗牙记录给出洗牙建议。

Agent 智能体	Role 职责
Multi-Expert Panel 专家会诊	3 role-prompted experts in PARALLEL (Pathologist on Claude; Epidemiologist + Dentist on Gemini) + deterministic Moderator. 三专家并行会诊 + 主持人汇总，保留分歧。
Consensus 共识核对	Second model re-reads the same image; agreement score. 第二个模型复核同一张图，计算一致性分数。
Retrieval (RAG) 知识检索	Grounds results in a curated oral-cancer knowledge base. 用口腔癌知识库为结论提供依据。
Patient Communication 患者沟通	Plain-language report with mandatory disclaimer. 通俗语言报告，强制附带免责声明。
Clinician Referral 转诊	Structured referral packet — generated ONLY on High risk. 仅在高风险时生成结构化转诊单。
Triage Prioritization 分诊排序	Urgency + SLA for the doctor queue. 为医生队列计算紧急度与时限。
Doctor Assistant 医生助手	REAL tool loop: model → tool → observation → repeat (≤4 steps) over a registry (profile / sessions / compare-all). 真实工具循环，最多 4 步，含全员比较工具。
Patient Chat 患者聊天	Multi-provider chain (Gemini → Claude → ChatGPT) + safe deterministic fallback. 多供应商链式回退 + 安全离线兜底。
Telegram Gateway 渠道网关	Zero-press channel adapter into the same pipeline and store. 零点击渠道适配器，与网页共用流水线与存储。

4. One Picture, Two Functions • 一图双检设计

A single shared prompt (lib/visionProviders/shared.ts) instructs whichever vision model is selected to return BOTH detections in one call: (1) oral-cancer cues — finding / region / probability / bounding boxes; (2) dental signs — cavity signs, gum redness, plaque/tartar, staining. The response is forced into typed JSON, coerced and clamped before any downstream agent reads it. The Dental Wellness agent then combines the image signs with toothbrush telemetry, the questionnaire, and the patient's last-scaling record to decide cavity risk and whether scaling is due.

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Safety: keys are never logged or echoed; the admin API returns masked previews only; the Gemini key travels in a header (not the URL) to keep it out of proxy logs.

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8. Safety & Limitations • 安全与限制

- Educational prototype only — not a medical device; every output carries a disclaimer; the chat agent refuses diagnosis requests.

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- Demo-grade auth (seeded plaintext-compared passwords) and demo-grade key storage (KV, not a secret manager) — replace both before any real use.

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- Risk weights are literature-calibrated heuristics, not a clinically validated model; the toothbrush telemetry is simulated pending real IoT hardware.

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- If a provider key is missing or a call fails, the pipeline degrades to mock and clearly labels the fallback in the provider badge and audit log.

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